



***Inami & Co., Ltd.***

**L-5110/5130**

**APPLANATION TONOMETER**





# INAMI APPLANATION TONOMETER L-5110/L-5130



**L-5110 INAMI Applanation Tonometer** for the measurement of intraocular pressure is designed to be mounted on the microscope permanently. When the tonometer is in use it can be swung forward just in front of the microscope and exact centering of the prism with the left eyepiece is engaged by notch. Illumination through the built-in cobalt filter projected to the plastic measuring prism for flattening the vertex of cornea gives reflection-free observation of split images.



**L-5130 INAMI Applanation Tonometer**, same as L-5110 Applanation Tonometer, can be used with various INAMI Slit Lamps.

L-5130 can be placed on the guide plate in one of two possible positions. These positions are connected with the microscope, for observation can be made either through the right or the left eyepiece.

## Specifications:

**Measuring range:** From 0 to 80mmHg in 2mmHg increments  $\pm 0.5\text{mmHg}$  **Diameter of the pneumatic face:** 3.06mm **Measurement:** 47mm wide  $\times$  30 mm deep  $\times$  85mm high **Weight:** 725gs.

## Standard Set

Applanation tonometer .....	1
Prism .....	2
Calibration bar .....	1
Special adaptor .....	1



Applanation Tonometer  
and prism set

*"A Passion for a Nice View"*

●Specifications subject to change without prior notice.



**Inami & Co., Ltd.**

No. 24-2, Hongo 3-Chome,  
Bunkyo-ku, Tokyo 113-0033, Japan  
Tel. 81-3-3814-5916 Fax. 81-3-5684-2126  
Website : <http://www.inami.co.jp>  
E-Mail : [corpinfo@inami.co.jp](mailto:corpinfo@inami.co.jp)

®TRADE MARK OF INAMI & CO.,LTD.



Certified Management System  
DIN EN ISO 9001:2000  
JIS Q 9001:2000  
EN ISO 13485:2000  
ISO 13485:2003



ZFG1000